

# DMC1 Antibody, HRP conjugated



PACO54479

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May participate in meiotic recombination, specifically in homologous strand assimilation, which is required for the resolution of meiotic double-strand breaks.

**Gene ID:**

DMC1

**Uniprot**

Q14565

**Synonyms:**

Meiotic recombination protein DMC1/LIM15 homolog, DMC1, DMC1H LIM15

**Immunogen:**

Recombinant Human Meiotic recombination protein DMC1/LIM15 homolog protein (5-153AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A