

# Pollen allergen Phl p 5a Antibody



PACO54506

---

## Product Information

**Size:**

50ug

**Reactivity:**

Phleum pratense

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

**Uniprot**

Q40962

**Synonyms:**

**Immunogen:**

Recombinant Phleum pratense Pollen allergen Phl p 5a protein (1-286AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A