

**Pollen allergen Phl p 5a Antibody, HRP
conjugated
PACO54507**



Product Information

Size:

50ug

Protein Background:

Reactivity:

Phleum pratense

Gene ID:

Source:

Rabbit

Uniprot

Q40962

Isotype:

IgG

Synonyms:

Applications:

ELISA

Immunogen:

Recombinant Phleum pratense Pollen allergen Phl p 5a protein (1-286AA).

Recommended dilutions:

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A