

ABCC2 Antibody, FITC conjugated



PACO54548

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Mediates hepatobiliary excretion of numerous organic anions. May function as a cellular cisplatin transporter.

Gene ID:

ABCC2

Uniprot

Q92887

Synonyms:

Canalicular multispecific organic anion transporter 1 (ATP-binding cassette sub-family C member 2) (Canalicular multidrug resistance protein) (Multidrug resistance-associated protein 2), ABCC2, CMOAT CMOAT1 CMRP MRP2

Immunogen:

Recombinant Human Canalicular multispecific organic anion transporter 1 protein (724-947AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A