

SLC30A3 Antibody



PACO54558

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Involved in accumulation of zinc in synaptic vesicles.

Gene ID:

SLC30A3

Uniprot

Q99726

Synonyms:

Zinc transporter 3 (ZnT-3) (Solute carrier family 30 member 3), SLC30A3, ZNT3


Immunogen:

Recombinant Human Zinc transporter 3 protein (286-388AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO54558 at dilution of 1:100.