

SLC30A3 Antibody, Biotin conjugated



PACO54561

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in accumulation of zinc in synaptic vesicles.

Gene ID:

SLC30A3

Uniprot

Q99726

Synonyms:

Zinc transporter 3 (ZnT-3) (Solute carrier family 30 member 3), SLC30A3, ZNT3

Immunogen:

Recombinant Human Zinc transporter 3 protein (286-388AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A