

ABCC3 Antibody, Biotin conjugated



PACO54613

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May act as an inducible transporter in the biliary and intestinal excretion of organic anions. Acts as an alternative route for the export of bile acid, and glucuronides from cholestatic hepatocytes.

Gene ID:

ABCC3

Uniprot

O15438

Synonyms:

Canalicular multispecific organic anion transporter 2 (ATP-binding cassette sub-family C member 3) (Multi-specific organic anion transporter D) (MOAT-D) (Multidrug resistance-associated protein 3), ABCC3, CMOAT2 MLP2 MRP3

Immunogen:

Recombinant Human Canalicular multispecific organic anion transporter 2 protein (827-953AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A