

CD3G Antibody, HRP conjugated



PACO54691

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

The CD3 complex mediates signal transduction.

Gene ID:

CD3G

Uniprot

P09693

Synonyms:

T-cell surface glycoprotein CD3 gamma chain (T-cell receptor T3 gamma chain) (CD antigen CD3g), CD3G, T3G

Immunogen:

Recombinant Human T-cell surface glycoprotein CD3 & gamma; chain protein (138-182AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A