

FGF17 Antibody



PACO54770

Product Information

Size:

50ug

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000,
IHC:1:20-1:200

Protein Background:

Plays an important role in the regulation of embryonic development and as signaling molecule in the induction and patterning of the embryonic brain. Required for normal brain development.

Gene ID:

FGF17

Uniprot

O60258

Synonyms:

Fibroblast growth factor 17 (FGF-17), FGF17


Immunogen:

Recombinant Human Fibroblast growth factor 17 protein (23-216AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Mouse liver tissue. All lanes: FGF17 antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 25, 24 kDa. Observed band size: 25 kDa.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO54770 at dilution of 1:100.