

FGF17 Antibody, Biotin conjugated



PACO54773

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays an important role in the regulation of embryonic development and as signaling molecule in the induction and patterning of the embryonic brain. Required for normal brain development.

Gene ID:

FGF17

Uniprot

O60258

Synonyms:

Fibroblast growth factor 17 (FGF-17), FGF17

Immunogen:

Recombinant Human Fibroblast growth factor 17 protein (23-216AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A