

GRK5 Antibody, FITC conjugated



PACO54808

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

GRK5

Uniprot

P34947

Synonyms:

G protein-coupled receptor kinase 5 (EC 2.7.11.16) (G protein-coupled receptor kinase GRK5), GRK5, GPRK5

Immunogen:

Recombinant Human G protein-coupled receptor kinase 5 protein (406-590AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A