

NR2F6 Antibody, FITC conjugated



PACO54912

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NR2F6

Uniprot

P10588

Synonyms:

Nuclear receptor subfamily 2 group F member 6 (V-erbA-related protein 2) (EAR-2), NR2F6, EAR2 ERBAL2

Immunogen:

Recombinant Human Nuclear receptor subfamily 2 group F member 6 protein (6-194AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A