

PTX3 Antibody, Biotin conjugated



PACO54969

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Plays a role in the regulation of innate resistance to pathogens, inflammatory reactions, possibly clearance of self-components and female fertility.

Gene ID:

PTX3

Uniprot

P26022

Synonyms:

Pentraxin-related protein PTX3 (Pentraxin-related protein PTX3) (Tumor necrosis factor alpha-induced protein 5) (TNF alpha-induced protein 5) (Tumor necrosis factor-inducible gene 14 protein) (TSG-14), PTX3, TNFAIP5 TSG14

Immunogen:

Recombinant Human Pentraxin-related protein PTX3 protein (99-164AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A