

SEPT3 Antibody, HRP conjugated



PACO54991

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Filament-forming cytoskeletal GTPase. May play a role in cytokinesis (Potential).

Gene ID:

37865

Uniprot

Q9UH03

Synonyms:

Neuronal-specific septin-3, SEPT3, SEP3

Immunogen:

Recombinant Human Neuronal-specific septin-3 protein (238-358AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A