

SLC25A12 Antibody, Biotin conjugated



PACO55013

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Catalyzes the calcium-dependent exchange of cytoplasmic glutamate with mitochondrial aspartate across the mitochondrial inner membrane. May have a function in the urea cycle.

Gene ID:

SLC25A12

Uniprot

O75746

Synonyms:

Calcium-binding mitochondrial carrier protein Aralar1 (Mitochondrial aspartate glutamate carrier 1) (Solute carrier family 25 member 12), SLC25A12, ARALAR1

Immunogen:

Recombinant Human Calcium-binding mitochondrial carrier protein Aralar1 protein (202-317AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A