

# TRIM3 Antibody



PACO55034

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Probably involved in vesicular trafficking via its association with the CART complex. The CART complex is necessary for efficient transferrin receptor recycling but not for EGFR degradation. Positively regulates motility of microtubule-dependent motor protein KIF21B.

**Gene ID:**

TRIM3

**Uniprot**

O75382

**Synonyms:**

Tripartite motif-containing protein 3 (Brain-expressed RING finger protein) (RING finger protein 22) (RING finger protein 97), TRIM3, BERP RNF22 RNF97

**Immunogen:**


Recombinant Human Tripartite motif-containing protein 3 protein (285-487AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO55034 at dilution of 1:100.