

VCP Antibody



PACO55054

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000,
IHC:1:20-1:200

Protein Background:

Gene ID:

VCP

Uniprot

P55072

Synonyms:

Transitional endoplasmic reticulum ATPase (TER ATPase) (EC 3.6.4.6) (15S Mg(2+)-ATPase p97 subunit) (Valosin-containing protein) (VCP), VCP


Immunogen:

Recombinant Human Transitional endoplasmic reticulum ATPase protein (704-806AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: Raw264.7 whole cell lysate. All lanes: VCP antibody at 3 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 90 kDa. Observed band size: 90 kDa.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO55054 at dilution of 1:100.