

HDC Antibody, FITC conjugated



PACO55068

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the biosynthesis of histamine from histidine.

Gene ID:

HDC

Uniprot

P19113

Synonyms:

Histidine decarboxylase (HDC) (EC 4.1.1.22), HDC

Immunogen:

Recombinant Human Histidine decarboxylase protein (504-648AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A