

BCAM Antibody, Biotin conjugated



PACO55085

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Laminin alpha-5 receptor. May mediate intracellular signaling.

Gene ID:

BCAM

Uniprot

P50895

Synonyms:

Basal cell adhesion molecule (Aubergier B antigen) (B-CAM cell surface glycoprotein) (F8/G253 antigen) (Lutheran antigen) (Lutheran blood group glycoprotein) (CD antigen CD239), BCAM, LU MSK19

Immunogen:

Recombinant Human Basal cell adhesion molecule protein (447-547AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A