

# USP10 Antibody, Biotin conjugated



PACO55097

---

## Product Information

|                               |   |
|-------------------------------|---|
| <b>Size:</b>                  | <b>Protein Background:</b>  |
| 50ug                          |   |
| <b>Reactivity:</b>            | <b>Gene ID:</b>   |
| Human                         | USP10   |
| <b>Source:</b>                | <b>Uniprot</b>  |
| Rabbit                        | Q14694  |
| <b>Isotype:</b>               | <b>Synonyms:</b>  |
| IgG                           | Ubiquitin carboxyl-terminal hydrolase 10 (EC 3.4.19.12) (Deubiquitinating enzyme 10) (Ubiquitin thioesterase 10) (Ubiquitin-specific-processing protease 10), USP10, KIAA0190 |
| <b>Applications:</b>          | <b>Immunogen:</b>   |
| ELISA                         | Recombinant Human Ubiquitin carboxyl-terminal hydrolase 10 protein (453-592AA).   |
| <b>Recommended dilutions:</b> | <b>Storage:</b>   |
|                               | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4  |

## Product Images

---

N/A

N/A