

SLC16A12 Antibody



PACO55186

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000,
IHC:1:20-1:200, IF:1:200-1:500

Protein Background:

Proton-linked monocarboxylate transporter that mediates creatine transport across the plasma membrane.

Gene ID:

SLC16A12

Uniprot

Q6ZSM3

Synonyms:

Monocarboxylate transporter 12 (MCT 12) (Creatine transporter 2) (CRT2) (Solute carrier family 16 member 12), SLC16A12, MCT12


Immunogen:

Recombinant Human Monocarboxylate transporter 12 protein (431-486AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4


Product Images



Western Blot. Positive WB detected in: Jurkat whole cell lysate, 293 whole cell lysate, HeLa whole cell lysate, HepG2 whole cell lysate. All lanes: SLC16A12 antibody at 3 μ g/ml. Secondary: Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 54 kDa. Observed band size: 54 kDa.



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using PACO55186 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO55186 at dilution of 1:100.