

ALYREF Antibody, HRP conjugated



PACO55207

Product Information

Size:	Protein Background:
50ug	
Reactivity:	Gene ID:
Human	ALYREF
Source:	Uniprot
Rabbit	Q86V81
Isotype:	Synonyms:
IgG	THO complex subunit 4 (Tho4) (Ally of AML-1 and LEF-1) (Aly/REF export factor) (Transcriptional coactivator Aly/REF) (bZIP-enhancing factor BEF), ALYREF, ALY BEF THOC4
Applications:	Immunogen:
ELISA	Recombinant Human THO complex subunit 4 protein (186-240AA).
Recommended dilutions:	Storage:
	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A