

IL17D Antibody



PACO55238

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Induces expression of IL6, CXCL8/IL8, and CSF2/GM-CSF from endothelial cells.

Gene ID:

IL17D

Uniprot

Q8TAD2

Synonyms:

Interleukin-17D (IL-17D) (Interleukin-27) (IL-27), IL17D


Immunogen:

Recombinant Human Interleukin-17D protein (16-202AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human gastric cancer using PACO55238 at dilution of 1:100.