

# LAMP3 Antibody, HRP conjugated



PACO55359

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

May play a role in dendritic cell function and in adaptive immunity.

**Gene ID:**

LAMP3

**Uniprot**

Q9UQV4

**Synonyms:**

Lysosome-associated membrane glycoprotein 3 (LAMP-3) (Lysosomal-associated membrane protein 3) (DC-lysosome-associated membrane glycoprotein) (DC LAMP) (Protein TSC403) (CD antigen CD208), LAMP3, DCLAMP TSC403

**Immunogen:**

Recombinant Human Lysosome-associated membrane glycoprotein 3 protein (28-224AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A