

# BOK Antibody, FITC conjugated



PACO55372

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

BOK

**Uniprot**

Q9UMX3

**Synonyms:**

Bcl-2-related ovarian killer protein (hBOK) (Bcl-2-like protein 9) (Bcl2-L-9), BOK, BCL2L9

**Immunogen:**

Recombinant Human Bcl-2-related ovarian killer protein (15-103AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A