

ABCC5 Antibody, FITC conjugated



PACO55380

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts as a multispecific organic anion pump which can transport nucleotide analogs.

Gene ID:

ABCC5

Uniprot

O15440

Synonyms:

Multidrug resistance-associated protein 5 (ATP-binding cassette sub-family C member 5) (Multi-specific organic anion transporter C) (MOAT-C) (SMRP) (pABC11), ABCC5, MRP5

Immunogen:

Recombinant Human Multidrug resistance-associated protein 5 protein (455-566AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A