

CA4 Antibody, FITC conjugated



PACO55416

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Reversible hydration of carbon dioxide. May stimulate the sodium/bicarbonate transporter activity of SLC4A4 that acts in pH homeostasis. It is essential for acid, overload removal from the retina and retina epithelium, and acid, release in the choriocapillaris in the choroid.

Gene ID:

CA4

Uniprot

P22748

Synonyms:

Carbonic anhydrase 4 (EC 4.2.1.1) (Carbonate dehydratase IV) (Carbonic anhydrase IV) (CA-IV), CA4

Immunogen:

Recombinant Human Carbonic anhydrase 4 protein (75-200AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A