

CPT2 Antibody, FITC conjugated



PACO55436

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

mitochondrial inner membrane, mitochondrion, nucleolus, nucleoplasm, carnitine O-palmitoyltransferase activity, carnitine shuttle, regulation of lipid metabolic process.

Gene ID:

CPT2

Uniprot

P23786

Synonyms:

Carnitine O-palmitoyltransferase 2, mitochondrial (EC 2.3.1.21) (Carnitine palmitoyltransferase II) (CPT II), CPT2, CPT1

Immunogen:

Recombinant Human Carnitine O-palmitoyltransferase 2, mitochondrial protein (388-551AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A