

FGFR4 Antibody, HRP conjugated



PACO55491

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

FGFR4

Uniprot

P22455

Synonyms:

Fibroblast growth factor receptor 4 (FGFR-4) (EC 2.7.10.1) (CD antigen CD334), FGFR4, JTK2 TKF

Immunogen:

Recombinant Human Fibroblast growth factor receptor 4 protein (25-182AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A