

FLOT1 Antibody, FITC conjugated



PACO55500

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Gene ID:

FLOT1

Uniprot

O75955

Synonyms:

Flotillin-1, FLOT1

Immunogen:

Recombinant Human Flotillin-1 protein (169-251AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A