

MST1 Antibody, HRP conjugated



PACO55651

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

extracellular exosome, extracellular region, extracellular space, receptor tyrosine kinase binding, hepatocyte growth factor receptor signaling pathway, negative regulation of gluconeogenesis, regulation of cAMP-dependent protein kinase activity.

Gene ID:

MST1

Uniprot

P26927

Synonyms:

Hepatocyte growth factor-like protein (Macrophage stimulatory protein) (Macrophage-stimulating protein) (MSP) [Cleaved into: Hepatocyte growth factor-like protein alpha chain; Hepatocyte growth factor-like protein beta chain], MST1, D3F15S2 DNF15S2 HGFL

Immunogen:

Recombinant Human Hepatocyte growth factor-like protein (363-465AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A