

NMBR Antibody



PACO55686

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:5000

Protein Background:

Receptor for neuromedin-B.

Gene ID:

NMBR

Uniprot

P28336

Synonyms:

Neuromedin-B receptor (NMB-R) (Epididymis tissue protein Li 185a) (Neuromedin-B-preferring bombesin receptor), NMBR


Immunogen:

Recombinant Human Neuromedin-B receptor protein (1-41AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Western Blot. Positive WB detected in: A549 whole cell lysate, HeLa whole cell lysate. All lanes: NMBR antibody at 4.4 μ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 44 kDa. Observed band size: 44 kDa.