

PNLIP Antibody, HRP conjugated



PACO55735

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

extracellular region, triglyceride lipase activity, lipid digestion, positive regulation of triglyceride lipase activity, retinoid metabolic process.

Gene ID:

PNLIP

Uniprot

P16233

Synonyms:

Pancreatic triacylglycerol lipase (PL) (PTL) (Pancreatic lipase) (EC 3.1.1.3), PNLIP

Immunogen:

Recombinant Human Pancreatic triacylglycerol lipase protein (60-214AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A