

SLC22A23 Antibody, HRP conjugated



PACO55787

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

integral component of membrane, organic anion transmembrane transporter activity.

Gene ID:

SLC22A23

Uniprot

A1A5C7

Synonyms:

Solute carrier family 22 member 23, SLC22A23, C6orf85

Immunogen:

Recombinant Human Solute carrier family 22 member 23 protein (3-206AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A