

# SLC22A23 Antibody, Biotin conjugated



PACO55789

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

integral component of membrane, organic anion transmembrane transporter activity.

**Gene ID:**

SLC22A23

**Uniprot**

A1A5C7

**Synonyms:**

Solute carrier family 22 member 23, SLC22A23, C6orf85

**Immunogen:**

Recombinant Human Solute carrier family 22 member 23 protein (3-206AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A