

# TNFAIP2 Antibody



PACO55806

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

May play a role as a mediator of inflammation and angiogenesis.

**Gene ID:**

TNFAIP2

**Uniprot**

Q03169

**Synonyms:**

Tumor necrosis factor alpha-induced protein 2 (TNF alpha-induced protein 2) (Primary response gene B94 protein), TNFAIP2

**Immunogen:**


Recombinant Human Tumor necrosis factor alpha-induced protein 2 protein (1-93AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in: Mouse liver tissue. All lanes: TNFAIP2 antibody at 6.9 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 73 kDa. Observed band size: 73 kDa.