

# POLD3 Antibody, Biotin conjugated



PACO55905

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

POLD3

**Uniprot**

Q15054

**Synonyms:**

DNA polymerase delta subunit 3 (DNA polymerase delta subunit C) (DNA polymerase delta subunit p66) (DNA polymerase delta subunit p68), POLD3, KIAA0039

**Immunogen:**

Recombinant Human DNA polymerase delta subunit 3 protein (203-275AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A