

CTF1 Antibody, HRP conjugated



PACO55923

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Induces cardiac myocyte hypertrophy in vitro. Binds to and activates the ILST/gp130 receptor.

Gene ID:

CTF1

Uniprot

Q16619

Synonyms:

Cardiotrophin-1 (CT-1), CTF1

Immunogen:

Recombinant Human Cardiotrophin-1 protein (1-201AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A