

SLC22A25 Antibody, Biotin conjugated



PACO55965

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

integral component of plasma membrane, inorganic anion exchanger activity, sodium-independent organic anion transmembrane transporter activity, urate transmembrane transporter activity, sodium-independent organic anion transport, urate transport.

Gene ID:

SLC22A25

Uniprot

Q6T423

Synonyms:

Solute carrier family 22 member 25 (Organic anion transporter UST6), SLC22A25, UST6

Immunogen:

Recombinant Human Solute carrier family 22 member 25 protein (31-145AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A