

# CD109 Antibody, FITC conjugated



PACO55968

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Modulates negatively TGFB1 signaling in keratinocytes.

**Gene ID:**

CD109

**Uniprot**

Q6YHK3

**Synonyms:**

CD109 antigen (150 kDa TGF-beta-1-binding protein) (C3 and PZP-like alpha-2-macroglobulin domain-containing protein 7) (Platelet-specific Gov antigen) (p180) (r150) (CD antigen CD109), CD109, CPAMD7

**Immunogen:**

Recombinant Human CD109 antigen protein (967-1122AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---

N/A

N/A