

CD109 Antibody, Biotin conjugated



PACO55969

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Modulates negatively TGFB1 signaling in keratinocytes.

Gene ID:

CD109

Uniprot

Q6YHK3

Synonyms:

CD109 antigen (150 kDa TGF-beta-1-binding protein) (C3 and PZP-like alpha-2-macroglobulin domain-containing protein 7) (Platelet-specific Gov antigen) (p180) (r150) (CD antigen CD109), CD109, CPAMD7

Immunogen:

Recombinant Human CD109 antigen protein (967-1122AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A