

BORA Antibody, Biotin conjugated



PACO55973

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Required for the activation of AURKA at the onset of mitosis.

Gene ID:

BORA

Uniprot

Q6PGQ7

Synonyms:

Protein aurora borealis (HsBora), BORA, C13orf34

Immunogen:

Recombinant Human Protein aurora borealis protein (178-312AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A