

COX15 Antibody, Biotin conjugated



PACO55981

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in the biosynthesis of heme A.

Gene ID:

COX15

Uniprot

Q7KZN9

Synonyms:

Cytochrome c oxidase assembly protein COX15 homolog, COX15

Immunogen:

Recombinant Human Cytochrome c oxidase assembly protein COX15 homolog protein (1-67AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A