

MACC1 Antibody, HRP conjugated



PACO55983

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Acts as a transcription activator for MET and as a key regulator of HGF-MET signaling. Promotes cell motility, proliferation and hepatocyte growth factor (HGF)-dependent scattering in vitro and tumor growth and metastasis in vivo.

Gene ID:

MACC1

Uniprot

Q6ZN28

Synonyms:

Metastasis-associated in colon cancer protein 1 (SH3 domain-containing protein 7a5), MACC1

Immunogen:

Recombinant Human Metastasis-associated in colon cancer protein 1 protein (371-514AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A