

# DDB2 Antibody



PACO56062

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human, Rat, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:500-1:5000

**Protein Background:**

**Gene ID:**

DDB2

**Uniprot**

Q92466

**Synonyms:**

DNA damage-binding protein 2 (DDB p48 subunit) (DDBb) (Damage-specific DNA-binding protein 2) (UV-damaged DNA-binding protein 2) (UV-DDB 2), DDB2

**Immunogen:**


Recombinant Human DNA damage-binding protein 2 protein (1-185AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

## Product Images

---



Western Blot. Positive WB detected in: Mouse liver tissue, Rat liver tissue. All lanes: DDB2 antibody at 4.9 $\mu$ g/ml. Secondary. Goat polyclonal to rabbit IgG at 1/50000 dilution. Predicted band size: 48, 27, 18, 41 kDa. Observed band size: 48 kDa.