

DUSP9 Antibody, FITC conjugated



PACO56124

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Inactivates MAP kinases. Has a specificity for the ERK family.

Gene ID:

DUSP9

Uniprot

Q99956

Synonyms:

Dual specificity protein phosphatase 9 (EC 3.1.3.16) (EC 3.1.3.48) (Mitogen-activated protein kinase phosphatase 4) (MAP kinase phosphatase 4) (MKP-4), DUSP9, MKP4

Immunogen:

Recombinant Human Dual specificity protein phosphatase 9 protein (142-263AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A