

MYOC Antibody, Biotin conjugated



PACO56129

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

MYOC

Uniprot

Q99972

Synonyms:

Myocilin (Myocilin 55 kDa subunit) (Trabecular meshwork-induced glucocorticoid response protein) [Cleaved into: Myocilin, N-terminal fragment (Myocilin 20 kDa N-terminal fragment); Myocilin, C-terminal fragment (Myocilin 35 kDa N-terminal fragment)], MYOC, GLC1A TIGR

Immunogen:

Recombinant Human Myocilin protein (183-294AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A