

OCLN Antibody, Biotin conjugated



PACO56353

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May play a role in the formation and regulation of the tight junction (TJ) paracellular permeability barrier. It is able to induce adhesion when expressed in cells lacking tight junctions.

Gene ID:

OCLN

Uniprot

Q16625

Synonyms:

Occludin, OCLN

Immunogen:

Recombinant Human Occludin protein (285-430AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A