

TBX3 Antibody, FITC conjugated



PACO56372

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Transcriptional repressor involved in developmental processes. Probably plays a role in limb pattern formation. Acts as a negative regulator of PML function in cellular senescence.

Gene ID:

TBX3

Uniprot

O15119

Synonyms:

T-box transcription factor TBX3 (T-box protein 3), TBX3

Immunogen:

Recombinant Human T-box transcription factor TBX3 protein (297-479AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images

N/A

N/A