

RASAL1 Antibody



PACO56458

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IF:1:50-1:200

Protein Background:

Probable inhibitory regulator of the Ras-cyclic AMP pathway. Plays a role in dendrite formation by melanocytes.

Gene ID:

RASAL1

Uniprot

O95294

Synonyms:

RasGAP-activating-like protein 1 (RAS protein activator like 1) (Ras GTPase-activating-like protein), RASAL1, RASAL


Immunogen:

Recombinant Human RasGAP-activating-like protein 1 protein (698-804AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, pH 7.4

Product Images



Immunofluorescence staining of MCF-7 cells with PACO56458 at 1:105, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).